

Sampling Technology	Collection Method	Blood Volume	Assay	Volumetric Plasma Collection	Storage and Transport	Suitability for Field Sampling	No Sample Prep Before Analysis	Non - Biohazard	Hematocrit Independent	No RBC Assay Interference
<u>Telimmune™</u>	Fingerstick or Lancet	25-60 µl	Plasma	Yes, ≈3 µl per Disc	Room Temp / Dry Mail	✓	✓	✓	✓	✓
DBS Cards	Fingerstick or Lancet	Typically ≥ 100 µl	Dried Whole Blood	X	Room Temp / Dry Mail	✓	X	✓	X	X
Whole Blood Microtube	Fingerstick or Lancet	≥ 5 µl	Plasma or Serum	X	≤ 42° F Liquid	Marginal	X Centrifuge	X	X	✓
Intravenous Blood Tube	Needles	≈ 20,000 µl	Plasma or Serum	Yes, but variable measurement	≤ 42° F Liquid	X	X Centrifuge	X	✓	✓

Green = Product Benefit, Red = Product Deficit, Yellow = In Between